

FlowConnector™ Kits - Researcher™

FlowConnector™ Kits - Researcher™ connect your WaterSep Researcher™ cartridge to your tubing and pump. The kits come with or without PendoTECH pressure sensors that enable monitoring of transmembrane pressure (TMP).




Diagram: FlowConnector™ Kit- Researcher™ without sensors

ORDERING INFORMATION

Place orders and request quotes online at www.watersep.net or email purchase orders to orders@watersep.net.

FlowConnector™ Kits - Researcher™ without sensors

FlowConnector™ Kit - Researcher™ Green

FlowConnector™ Kit - Researcher™ Steamer

FlowConnector™ Kit - Researcher™ ReUse

Product Code: FC RES KIT-000 01

Includes:

For Retentate/Feed

Item 1: (2) 3/4" TC Clamps

Item 2: (2) 3/4" TC x 3/8" Barb (Reducers)

Item 3: (4) 3/4" TC Silicone Gaskets (2 extra not shown)

For Retentate

Item 5: (1) 3/4" Pinch Valve

For Permeate

Item 4: (1) 3/8" Pinch Valve

Item 6: (2) Press In Plugs, 3/8" Barb (1 extra not shown)

FlowConnector™ Kits - Researcher™




Diagram: Flow Connector Kit- Researcher™ with sensors

FlowConnector™ Kits - Researcher™ with sensors

FlowConnector™ Kit - Researcher™ Green

FlowConnector™ Kit - Researcher™ Steamer

FlowConnector™ Kit - Researcher™ ReUse

Product Code: FC RES KIT-SEN 01

Includes:

For Retentate/Feed

Item 1: (2) 3/4" TC Clamps

Item 2: (4) 3/4" TC Silicone Gaskets (2 extra not shown)

Item 7: (2) Single Use Pressure Sensor 3/4" TC x 3/8" Barb

For Retentate

Item 4: (1) 3/4" Pinch Valve

For Permeate

Item 3: (1) 3/8" Pinch Valve

Item 5: (2) Press In Plugs, 3/8" Barb (1 extra not shown)

Item 6: (1) Single Use Pressure Sensor 1/4" Barb

MATERIALS OF CONSTRUCTION FOR KIT COMPONENTS

IN FLUID CONTACT

Gaskets:

Platinum-cured silicone

Reducer and Plugs:

Polypropylene AF

Sensors:

Polysulfone

NOT IN FLUID CONTACT

Pinch Valves:

Glass Reinforced White Nylon

Clamps:

Glass Reinforced White Nylon

Technical Support: experttalk@watersep.net
or call 508.970.0089 x 204